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Dermotrace® **Instructions for Use**

Dermotrace® is a sterile skin marker for surgical use.

Presentation

Three types with different tips are available.

- Dermotrace® bi-pointe with two tips: one medium tip (4.5 mm diameter) for a line thickness of 2 to 4 mm and one fine tip (2 mm diameter) for a line thickness of 1 to 2 mm.
- Dermotrace® monopointe fine has one fine tip (2 mm diameter) for a line thickness of 1 to 2 mm.
- Dermotrace® monopointe extra-fine has one extra-fine tip (0.7 mm diameter) for a line thickness of 0.4 to 0.9 mm.

Facilitating surgeons' work, the Dermotrace® skin marker comes sterile in a set with a flexible graduated ruler of 15 cm length.

Precautions for use

- Apply the marker to clean and dry skin only, preferably before disinfection. In case the line drawn before disintection is blurred by the antiseptic, redraw using a new Dermotrace® marker. If the marking cannot be made before disinfection, it is recommended to use Dermotrace® bi-pointe or Dermotrace® monopointe fine, as their ink capacity is larger than that of Dermotrace® monopointe extra-fine.
- On very dark skin, Dermotrace® is not well visible.
- Dermotrace® may not be used with patients having an allergy against one of the components.
- The sterile pouch should be opened above a sterile field, avoiding the risk of Dermotrace® falling onto a contaminated surface.

Single Use

Dermotrace® is for single use only.

There is a risk of contamination of skin markers with pathogenic germs and / or commensal germs. Therefore it is recommended to use several markers for marking different parts of the body.

Technical features

- Perfect ink resistance to soap and water
- Good visibility with most skin types
- Can be erased easily with alcoholic solution
- Easy drawing even in angles and wrinkles thanks to the fine tip
- Medical grade ink, biocompatible and nontoxic

Ink composition

Methyl violet, ricinus oil, diethylen glycol, osmotic water

Disposal

Disposal according to national regulations for waste and / or - after patient contact or contact with contaminated material - to regulations for contaminated waste.

REF	Description Sterile marker, together with a ruler in a peelable pouch	Items per dispenser box
18996	Dermotrace® bi-pointe: medium tip, Ø 4.5 mm and fine tip, Ø 2 mm	50
18986	Dermotrace® monopointe, extra-fine tip, Ø 0.7 mm	50
18976	Dermotrace® monopointe, fine tip, Ø 2 mm, transparent cap	50

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